

Bachelor Thesis 2016

Radiosensitization using gold nanoparticles


Lies Deceuninck en Hannelore Verhoeven

Assistents: Bert De Roo
Mattias Vervaele
Professor: Chris Van Haesendonck

DNA damage using ionizing radiation

- Chemotherapy
- Surgery
- **Radiation therapy**


Energy \sim MeV




Source: <https://chem.kuleuven.be/veiligheid/info/ioniserende-st.htm>

Radiosensitization of cancer cells with gold nanoparticles (GNP) $E \sim \text{keV}$

Photoelectric absorption




Compton effect



Targeting of the GNP to the tumor


[c] **Passive targeting**

PEG coating







Active targeting

Antibodies



1. Synthesis
2. Characterization
3. Radiosensitization




Reduction of gold ions to form GNP




Gold ions: HAuCl_4 solution

Reducing agent: $\text{Na}_3\text{C}_6\text{H}_5\text{O}_7$

Reduction of gold ions to form GNP






The amount of citrate controls the size




Citrate 1%

100ml HAuCl₄ 0.01%

PEG for targeting and stabilization




20k, 10k, 5k, 1k



TEM image analysis to determine core diameter

UV-Vis spectroscopy


1. Add PEG
2. Size GNP
3. Add NaCl
4. Size GNP



bigger size → too little PEG
same size → enough PEG


Results

GNP no PEG




Results

15nm GNP 20k PEG for different PEG/GNP




Results

15nm GNP 20k PEG for different PEG/GNP



Results

15nm GNP 20k PEG for different PEG/GNP



Overview

Introduction

Synthesis GNP

Chemical Protocol

Size GNP

Stabilization

Characterization

Size GNP


Chemical Protocol

UV-VIS


TEM

Hydrodynamic Radius



DLS



Dynamic light scattering (DLS)




Hydrodynamic radius (R_h)
 → Rayleigh scattering




Results

Functionalisation no PEG



Results

Functionalisation 20k PEG




Results

Functionalization 15nm 20k PEG

Proportion (PEG/GNP)	Average
5/10	51.93 ± 2.76
6/10	80.89 ± 14.64
7/10	65.24 ± 14.32
8/10	83.91 ± 18.42
9/10	

Original functionalization 20k (8/10)




Results

Functionalization 15nm 20k PEG

Proportion (PEG/GNP)	Average	Average (centrifuge)
5/10	51.93 ± 2.76	68.70 ± 7.99
6/10	80.89 ± 14.64	65.16 ± 11.61
7/10	65.24 ± 14.32	57.73 ± 7.72
8/10	83.91 ± 18.42	72.36 ± 10.44
9/10		56.54 ± 3.91

Original functionalization 20k (8/10)



Conclusion

- Synthesis of GNP
- Characterization
- Stabilization with neutral PEG
- Stabilization with positively charged PEG
- X-Rays
- Analyze effect on DNA
- Solve problem with DLS

