




# WALLIS - Zeta Potential Analyzer




*An introduction to electrophoresis principle  
for particle charge measurements*






# Why measuring Zeta potential?



***Charged particles***  
repel each other




***Uncharged particles***  
are free to collide and  
aggregate


## Colloidal solution stability ?



Stable systems



Flocculations



Coagulation




Sedimentation

- ➔ **Electrostatic or charge stabilization** : This effect uses more natural interactions between particles through the distribution of charged species (ions) in the solution.





# Zeta Potential definition




- If surface has - charge, then Cations are attracted to surface
- Anions attracted to Cations, builds electric double layer
- Slipping plane: distance from particle surface where ions move with particle
- Zeta Potential = potential (mV) at slipping plane





# Electrophoresis Effect



*Applying an Electric field: Only strongly attached ions move with the particle.*

*Particle velocity relies directly to its electrophoretic mobility under a known E-field*






# Laser Doppler Electrophoresis (LDE) principle

The Idea:

- measure the speed of the charged nano particles by applying a known electric fields
- use a laser and the Doppler shift to measure the averaged speed of the particles


## Heterodyne optical Interferometer





# Measurement Sequence

## Heterodyning Signal



*FFT*

Doppler shift

$$f_D$$

Measured

Electrophoretic  
Mobility

$$\mu_e$$

By *LDE*  
(Laser Doppler  
Electrophoresis)

Double Layer model

Huckel ?  
Smoluchowski ?

$f(\kappa \cdot a)$   
Henry function

Computed

Zeta  
Potential

$$\zeta$$





# Measurement of mobility

## Electrophoretic mobility

$$\mu_e = \frac{\lambda}{E \sin \theta_d} f_D$$


Where :

$\mu_e$ : electrophoretic mobility ( $\mu\text{mcm/Vs}$ )

$E$  : applied electric field

$f_{las}$  : laser frequency

$f_D$  : Doppler frequency



Particle motion causes Doppler shift

Frequency -> mobility





# Zeta Potential : Calculation of ZETA

## Zeta potential

$$\zeta = \mu_e \frac{\eta}{\varepsilon} f(\kappa \cdot a)$$

Where :

$\zeta$ : Zeta Potential (mV)

$\mu_e$ : electrophoretic mobility

$\varepsilon$ : medium permittivity


$\eta$ : medium viscosity

$a$ : particle radius

$\kappa$ : invert of double layer thickness

$\kappa^{-1}$ : Debye length

$f(\kappa \cdot a)$ : Henry's function



Zeta is dependent of the  $\kappa \cdot a$  factor value

$\kappa$  is dependent of the solvent








# Zeta Potential : Fundamental equation

General cases



$$1 \leq f(\kappa \cdot a) \leq 1.5$$


Aqueous solution



$$f(\kappa a) = 1$$

Smoluchowski

Organic solution



$$f(\kappa a) = 1,5$$

Hückel

99% cases !!!

Calculated

Measured


$$\zeta = \mu_e \frac{\eta}{\varepsilon} f(\kappa \cdot a)$$

Needs to know the value !

Zeta potential calculation is based on the knowledge of the double layer, define by Henry function  $f(\kappa \cdot a)$ .



# *WALLIS $\zeta$ : High-Resolution Zeta Potential Analysis*



- Charge/zeta potential measurement of nano-particles in suspension
- Complementary tool to VASCO for colloid characterization
- High resolution analysis from -200mV to + 200mV






# Direct acquisition : No correlator for higher resolution

With correlator  
(other supplier)

correlator

Averaging




FFT




Maximum time  
*unlimited*


FFT



Averaging



Without High  
resolution acquisition  
data board (Wallis)






# Electrode Assembly


Capillary = 3 mm

Electro-osmosis effect



Electrode

Charged interface



Dip cell = 10 mm

Connector

Hellma Cell

Quartz, Glass or  
polystyrene polished

Vitreous Carbon  
Electrodes

Laser beam

Min. sample  
volume  
750  $\mu$ L

The Dip cell small electrodes configuration is based on established principles.  
No electro-osmosis in WALLIS cell ▶ The need to focus at any "stationary plane" is eliminated.





# Electrode Assembly

Concept oriented for handling easiness



- Easy filling from above !
  - No risk of bubbles ! (vs capillary cell)
- 
- Electrodes are independent (moderate cost)
  - Vitreous carbon electrodes (lifetime > 1 year - extreme resistance to chemical attack - low electrical resistance )
  - Conventional cleaning of electrodes (ultra sound or bath tub)
- 
- Life of a quartz cell (> 5 years)
  - Excellent optical quality (quartz and / or polystyrene polished !)
  - Compatible with organic solvent
  - Compatible with Hellma cell ...






## Example : sample of industrial surfactant

High resolution of Wallis allow high quality measurement for applications:

- Precise measurement of IEP
- Functionalization of colloid
- Repeatability of particles synthesis
- Sharp variation in particles stability










# ZetaQ

- Unique and proprietary software
- User friendly and touch-screen oriented
- Database structure





# Thank you of your attention!

**Customer support:**

**Boris PEDRONO**

**Mail:** [boris.pedrono@cordouan-tech.com](mailto:boris.pedrono@cordouan-tech.com)

**CORDOUAN TECHNOLOGIES**

**+33 556 157 541**

**Cité de la photonique**

**11, Avenue Canteranne**

**33600 Pessac FRANCE**

**[www.cordouan-tech.com](http://www.cordouan-tech.com)**

**[Follow Cordouan Technologies on Twitter!](#)**

