


Bachelor Thesis 2016

Radiosensitization using gold nanoparticles

Lies Deceuninck en Hannelore Verhoeven

Assistents: Bert De Roo
Mattias Vervaele
Professor: Chris Van Haesendonck




- Chemotherapy
- Surgery
- **Radiation therapy**

The diagram illustrates two pathways of DNA damage. On the left, a DNA double helix is shown with base pairs (G-C, A-T, T-A, G-C, G-C, T-A, T-A, A-T, A-T, T-A, C-G, G-C, T-A, G-C, T-A, C-G, T-A, A-T, C-G, T-A) and arrows indicating the direction of the strands. Two black dots represent sites of damage on the DNA. To the right, the 'Direct' pathway is shown, where a wavy line representing ionizing radiation hits a DNA molecule, with an inset showing an electron (e-) and a proton (p+) pair. Below this, the 'Indirect' pathway is shown, where H₂O is converted to HO (hydroxyl radical) by radiation, with an inset showing an electron (e-) and a proton (p+) pair.


Katholieke Universiteit Leuven

Radiosensitization of cancer cells with gold nanoparticles (GNP) $E \sim \text{keV}$

Photoelectric absorption




Compton effect



Targeting of the GNP to the tumor


Passive targeting

PEG coating



Active targeting





Antibodies




Overview Project

Radiosensitization of cancer cells using gold nanoparticles

1. Synthesis
2. Characterization
3. Radiosensitization





Reduction of gold ions to form GNP




Gold ions: HAuCl_4 solution

Reducing agent: $\text{Na}_3\text{C}_6\text{H}_5\text{O}_7$

Reduction of gold ions to form GNP






The amount of citrate controls the size




Citrate 1%

100ml HAuCl₄ 0.01%



PEG for targeting and stabilization




20k, 10k, 5k, 1k



TEM image analysis to determine the core diameter




TEM image analysis to determine the core diameter



Exp. Size (nm)	Size (nm)
15	12.98 ± 0.23
	2.99 ± 0.16
30	18.29 ± 0.23
45	46.75 ± 0.47

UV-Vis spectroscopy


1. Add PEG
2. Size GNP
3. Add NaCl
4. Size GNP



bigger size \rightarrow too little PEG
same size \rightarrow enough PEG


Results

GNP no PEG




Results

15nm GNP 20k PEG for different PEG/GNP




Results

15nm GNP 20k PEG for different PEG/GNP



Results

15nm GNP 20k PEG for different PEG/GNP



Overview

Introduction

Synthesis GNP

Chemical Protocol

Size GNP

Stabilization

Characterization

Size GNP


Chemical Protocol

UV-VIS


TEM

Hydrodynamic Radius



DLS



Dynamic light scattering (DLS)




Hydrodynamic radius (R_h)
 \rightarrow Rayleigh scattering




Results

Functionalisation no PEG



Results

Functionalisation 20k PEG




Results

Functionalization 15nm 20k PEG

Proportion (PEG/GNP)	Average
5/10	51.93 ± 2.76
6/10	80.89 ± 14.64
7/10	65.24 ± 14.32
8/10	83.91 ± 18.42
9/10	

Original functionalization 20k (8/10)




Results

Functionalization 15nm 20k PEG

Proportion (PEG/GNP)	Average	Average (centrifuge)
5/10	51.93 ± 2.76	68.70 ± 7.99
6/10	80.89 ± 14.64	65.16 ± 11.61
7/10	65.24 ± 14.32	57.73 ± 7.72
8/10	83.91 ± 18.42	72.36 ± 10.44
9/10		56.54 ± 3.91

Original functionalization 20k (8/10)



Conclusion

- Synthesis of GNP
- Characterization
- Stabilization with neutral PEG
- Stabilization with positively charged PEG
- X-Rays
- Analyze effect on DNA
- Solve problem with DLS

